

THE LITTLE RICHIE BRACE®

An alternative to traditional AFO therapy for the pediatric patient



Modeled after the successful Richie Brace® for adults, the Little Richie Brace® combines a custom balanced functional foot orthotic with a semi-rigid ankle brace.

FEATURES

UCBL style custom balanced
Foot orthosis
35 mm heel cup, medial and lateral
flanges

Custom, individualized placement
Of anatomic ankle hinge

Custom semi-rigid limb uprights
On sides of lower leg

Aesthetic colors and design

Full flexion and fixed ankle pivot
Options

BENEFITS

Intrinsic control of midtarsal
and subtalar joints
Transverse and Frontal plane
control of the foot

Pure sagittal plane ankle
movement

Low bulk, sports style brace

Non-medical sports appearance

Versatile clinical applications

CLINICAL INDICATIONS

for children age 1 to age 12

Neuromuscular Conditions: Cerebral palsy, Muscular Dystrophy,
Ataxia, Hypotonia

Congenital Deformities: Talipes equino varus, Talipes calcaneal
Valgus, Tarsal coalition

Casting is Required—Please follow Richie Brace® casting instructions* and provide a mid-calf measurement of circumference. Or, use STS AFO casting sock***

Suggested Reimbursement Billing Codes: STS Casting sock **A4590**
AFO, molded to patient model with Ankle joint **L1970** (use modifier RT or LT)
Soft interface, below knee **L2820**