

## RICHIE BRACE ® TREATMENT GUIDE

					
<b>Richie Brace® Standard</b>	<b>Dynamic Assist</b>	<b>Restricted Hinge Pivot</b>	<b>Restricted Hinge w/Arch Suspender</b>	<b>California AFO</b>	<b>Richie Solid AFO</b>
Stage II PTTD	Dropfoot w/out equinus No spasticity Stable knee  (must have all 3 above)	Dropfoot w/ equinus  Dropfoot w/spasticity  Stable knee	Medial: Stage II or III PTTD with subluxed T-N Joint	Severe Deformity: Stage IV PTTD  Charcot Deformity	Dropfoot with unstable knee
Chronic Ankle Instability	Post CVA  Peroneal nerve injury  Post Polio	DJD of Hindfoot and Ankle	Lateral: Peroneal tendinopathy Fixed varus deformity of hindfoot/ankle	Severe DJD of Ankle or Hindfoot	
<b>C A S T I N G</b>					
Ankle Casting Sock	Ankle Casting Sock	Ankle Casting Sock	Ankle Casting Sock	Mid Leg Casting Sock	Full Leg (Bermuda) Casting Sock