

INTRODUCING

THE *Richie* GAUNTLET™

Legitimate Varus/Valgus control of the hindfoot



RICHIE GAUNTLET WITH LATERAL ARCH SUSPENDER



RICHIE GAUNTLET WITH MEDIAL ARCH SUSPENDER

FEATURES:

Varus or Valgus lifting strap • Intrinsic balanced AFO shell • Extrinsic rearfoot post

INDICATIONS:

Rigid, non-reducible Adult Acquired Flatfoot (Stage III & IV) • Severe DJD or Deformity of hindfoot • Charcot Arthropathy

CASTING INSTRUCTIONS:

Cast for the Richie Gauntlet using the STS Mid Leg Sock



Non-Weight & Weight bearing casts are accepted

Available in 2 colors, Tan & Chocolate



2 shell heights: 7" (left) and 9" (right)

BILLING CODES:

(* Ceiling & Floor Approx. Fees, Region D)

L1940 AFO, Custom-fabricated (*Region D \$517 - \$388)

L2275 Varus/Valgus Control (*Region D \$131 - \$98)

L2280 Molded inner boot (*Region D \$474 - \$355)

L2820 Soft interface (*Region D \$90 - \$68)

AVAILABLE THROUGH:

BURNS INTERNATIONAL, INC.

Orthotic Laboratories

Facilities in Nebraska and New Jersey

1-800-262-8767